

Sygnature Discovery Ltd, BioCity, Pennyfoot Street, Nottingham, NG1 1GF, United Kingdom.

The new libraries being synthesised will focus on: **novelty**, **shape** (eg. sp³-rich cores), **diversity potential** and **innovative library design** and will complement classical compound collections.

- cLogP <3
- MW <300
- 3D shape required
- Preferably amine, acid or alcohol handle
- No reactive groups in final compounds (e.g. imine, Michael acceptor...)



A selection of cores generated at
Sygnature Discovery

- Solubility
- Stability
- Ionisation
- Absorbency
- Reactivity
- Isolation
- Purification
- Recovery

Amide coupling

$\text{HNR}_2 + \text{R}^1\text{COOH} \xrightarrow[\text{DMF}]{\text{HATU, DIPEA}} \text{R}^1\text{CON}(\text{R})_2$

Reductive amination

$\text{HNR}_2 + \text{R}^1\text{CHO} \xrightarrow[\text{DMF}]{\text{STAB, AcOH}} \text{R}^1\text{CH}_2\text{N}(\text{R})_2$

Sulfonylation

$\text{HNR}_2 + \text{R}^1\text{SO}_2\text{Cl} \xrightarrow[\text{DMF}]{\text{pyridine}} \text{R}^1\text{SO}_2\text{N}(\text{R})_2$

Boc deprotection

$\text{R}_2\text{N-CO-O-C(CH}_3)_3 \xrightarrow[\text{DCM}]{\text{TFA}} \text{R}_2\text{N-H}$

Example library reactions tested during validation

Production and Dispatch

```

graph LR
    SM[Starting materials  
8-10 Scaffolds per idea  
Reagent sets are stored as stock solutions in DMA in barcoded 48 well plates] -- Plate up  
Reaction mixtures prepared using liquid handler or manual pipette --> RP[Reaction plate  
• 96 Well reaction plate  
• 96 Reactions per scaffold  
• 0.1 mmol scale  
• ~800 Reactions total  
• Shake at ambient temperature overnight]
    RP --> FP[Filter plate  
Reaction mixtures filtered prior to purification]
    FP -- Purification  
Semi-prep HPLC can purify >100 samples a day --> BT1[Barcoded tubes  
Fractions collected by semi-prep HPLC into barcoded tubes and dried by centrifugal evaporation]
    BT1 -- Reformat  
Liquid handler --> AP1[Analysis plate]
    AP1 -- UPLC analysis, weighing --> SV1[Submission vials  
Dry compounds re-dissolved in MeOH/CHCl3 on the liquid handler and transferred to vials for ELF and for Sygnature Discovery  
Compounds dried for a final time]
    SV1 -- Dispatch  
Compounds passing QC are delivered to Bioascent, Scotland --> CH[Compound hub  
Compounds made into 10 mM solutions in DMSO and added to screening collection]
    AP1 -- UPLC analysis, weighing --> SV2[Submission vials  
QC pass criteria  
Purity >85% by UPLC  
Quantity 5-50 μmol]
    SV2 --> CH
    
```

Reagent set

Waters UPLC

Analysis plate

Analysis plate
Samples diluted with DMSO into an analysis plate using the liquid handler

UPLC analysis
Plate stacker enables us to queue up to 9 plates

Analysis results
Focussed gradients for purification automatically generated from UPLC retention times

Filter plate
Reaction mixtures filtered prior to purification

Purification
Semi-prep HPLC can purify >100 samples a day

Barcoded tubes
Fractions collected by semi-prep HPLC into barcoded tubes and dried by centrifugal evaporation

Reformat
Liquid handler

Analysis plate

UPLC analysis, weighing

Submission vials
Dry compounds re-dissolved in MeOH/CHCl₃ on the liquid handler and transferred to vials for ELF and for Sygnature Discovery
Compounds dried for a final time

QC pass criteria
Purity >85% by UPLC
Quantity 5-50 μmol

Dispatch
Compounds passing QC are delivered to Bioascent, Scotland

Compound hub
Compounds made into 10 mM solutions in DMSO and added to screening collection

GeneVac HT6 centrifugal evaporator

Starting materials
8-10 Scaffolds per idea
Reagent sets are stored as stock solutions in DMA in barcoded 48 well plates

Plate up
Reaction mixtures prepared using liquid handler or manual pipette

Reaction plate

- 96 Well reaction plate
- 96 Reactions per scaffold
- 0.1 mmol scale
- ~800 Reactions total
- Shake at ambient temperature overnight

Plating up using reagent sets

Reaction plate

Waters semi-preparative HPLC

Tecan liquid handler

MW v cLogP of scaffolds for first 20,000 compounds

Fsp³ distribution of scaffolds