

Introduction

Transient Receptor Potential (TRP) channels constitute a vital group of integral membrane proteins that are widely distributed in mammalian cells and tissues. These channels play a crucial role in defining cellular signalling and function through regulation of membrane excitability and mediation of intracellular calcium levels. While TRP channels present promising therapeutic opportunities, the identification of compounds devoid of cross-reactivity with closely related family members remains a significant challenge.

TRPML1 and its Multifaceted Roles

TRPML1, a non-selective cation-permeable channel, is predominantly found on the membranes of late endosomes and lysosomes across mammalian cell types. TRPML1 is pivotal in maintaining lysosomal calcium homeostasis, modulating autophagy, and influencing oxidative stress. Mutations in TRPML1 have been implicated in various conditions, including lysosomal storage disorders and neurodegenerative diseases such as Alzheimer's and Parkinson's disease. Correcting aberrant TRPML1 activity has been suggested as a potential strategy to rectify abnormal autophagy and help clear protein/ROS accumulation in affected cells.

TRP Channel Discovery Platform

Despite some success in modulating TRPML1 with synthetic compounds, the lack of specific pharmacological tools has hampered a deeper understanding of TRPML1's role in normal lysosomal function and its association with pathological processes.

SB Drug Discovery's comprehensive TRP channel discovery platform includes meticulously optimized TRP channel cell lines, including TRPML1, TRPML2, and TRPML3, coupled with robust assays designed to support your drug discovery research. Leveraging both 384-well fluorescence-based and electrophysiology platforms for high-throughput screening and lead optimization, SB Drug Discovery enables rapid hit identification and downstream discovery cascades.

TRP Channel Discovery Platform

•Comprehensive TRP Channel Expertise

Broad scientific expertise, optimized TRP channel cell lines and assays, facilitating high-throughput screening and lead optimization.

•Robust optimized TRPML1 Assays

Tightly controlled assays to ensure reproducibility and reliability, aiding hit finding and downstream electrophysiology studies.

•Rapid Design-Make-Test-Analyze Support

SB Drug Discovery's TRP channel platform supports rapid cycles of design-make-test-analyze to expedite TRP channel drug discovery.

Methods

HEK cells stably expressing TRPML1 were produced by SB Drug Discovery. Whole cell patch-clamp experiments were carried out at room temperature using multi-hole chips on the SynchroPatch 384i automated electrophysiology platform. Currents were elicited by using repeated ramps, steadily increasing from -140 mV to +100 mV over 190 ms, from a holding potential of 0 mV. Data analysis was performed using Data Control 384 V2.3 (Nanion) and GraphPad Prism V10.1.





Figure 1. TRPML1 Voltage protocol.

TRPML1 currents were elicited by using repeated ramps, steadily increasing from -100 mV to +100 mV over 190 ms, from a holding potential of 0 mV.

TRPML1 Electrophysiology

To validate the TRPML1 electrophysiology assay, reference agonist ML SA5 was assessed at a range of concentrations using the SyncroPatch384 automated electrophysiology platform. Figure 2A highlights the success rate of the assay, with > 85% of cells successfully completing the experiment and passing strict QC criteria. The superimposed current traces and corresponding time courses shown in Figures 2B and 2C show modulation of TRPML1 activity in the presence of increasing concentrations of ML SA5 (0.01-10 μ M). The EC₅₀ of ML SA5 was 0.5 μ M which is in line with literature.



TRPML1 Pharmacology

SB Drug Discovery's TRPML1 assays have undergone extensive characterization and can be conducted in agonist, antagonist, or dual mode, providing a comprehensive understanding of compound effect and potency. Each assay meets the highest standard of quality control by integrating standard reference compounds, generating robust and reliable data.





Figure 3. TRPML1 pharmacology.

Concentration response curves comparing two reference agonists (A), ML SA1 and ML SA5 and three reference antagonists (B), GW 405833, Lanthanum Chloride and ML SI3 against the TRPML1 channel. Data is generated using the automated electrophysiology platform.

Platform correlation

Fluorescence-based high-throughput screening is the most commonly used route for ion channel hit identification. SB Drug Discovery's TRPML1 high-throughput screening assay displays excellent correlation with downstream automated electrophysiology assays, ensuring seamless transition from hit identification to lead optimization.



Reference Compound	Platform	EC ₅₀ (μM)
ML SA1	Fluorescence	4.8
	Electrophysiology	5.6
ML SA5	Fluorescence	0.6
	Electrophysiology	0.6
GW405833	Fluorescence	1.3
	Electrophysiology	1.9

Figure 4. Platform correlation.

Concentration response curves for reference compounds (A) ML SA1, (B) ML SA5, and (C) GW 405833 generated using fluorescence and automated electrophysiology platforms displaying excellent cross-platform correlation.

TRPML1 Screening Cascade

Enabling rapid design-make-test-analyze cycles, SB Drug Discovery's TRP Channel screening cascade helps streamline the identification, confirmation and optimization of potential drug candidates, offering a fast-track solutions in the pursuit of novel therapeutics.




Figure 5. SB Drug Discovery's TRPML1 screening cascade.

TRP Channel Selectivity Panel

At SB Drug Discovery, we recognize the critical importance of accessing a diverse array of TRP channel assays to support your discovery research. Our TRP channel discovery platform boasts a comprehensive suite of highly optimized assays enabling rapid screening of libraries and in-depth profiling of your compounds for TRP channel cross-reactivity, allowing you to identify compounds of interest and assess selectivity throughout your discovery campaign.

Human TRP Channel Collection				
TRPA1	TRPC1	TRPM2	TRPML1	TRPV1
	TRPC3	TRPM3	TRPML2	TRPV2
	TRPC4	TRPM4	TRPML3	TRPV3
	TRPC5	TRPM5		TRPV4
	TRPC6	TRPM8		TRPV5
	TRPC7			TRPV6

TRP Channel Species Variants			
TRPA1 (mouse)	TRPC5 (mouse)	TRPV1 (mouse)	TRPV4 (mouse)
TRPA1 (rat)	TRPC5 (rat)	TRPV1 (rat)	TRPV4 (rat)
TRPA1 (guinea pig)	TRPC6 (mouse)	TRPV2 (mouse)	TRPV5 (mouse)
TRPA1 (mini pig)	TRPC6 (rat)	TRPV2 (rat)	TRPV5 (rat)
TRPA1 (dog)	TRPM3 (rat)	TRPV3 (mouse)	TRPV6 (mouse)
TRPA1 (sheep)		TRPV3 (rat)	TRPV6 (rat)
TRPA1 (monkey)			

SB Drug Discovery's Commitment to Ion Channel Research

SB Drug Discovery scientists leverage their extensive experience in ion channel drug discovery to expedite the identification of novel TRPML1 modulators through rigorous high-throughput screening assays and premium-quality reagents. Our meticulously validated TRPML1 screening assays are ready for immediate implementation, allowing for fast turnaround times and the delivery of precise, reproducible data to advance your drug discovery research.



Cell Line Generation



Drug Discovery Screening



Ion Channel Electrophysiology



Inflammation

SB Drug Discovery
Telford Pavilion, Todd Campus
West of Scotland Science Park
Glasgow
G20 0XA
United Kingdom

T: +44 (0)141 587 6100
enquiries@sdrugdiscovery.com
sdrugdiscovery.com

